

## ACTIVITY OF PERA SAFE™ AGAINST *BACILLUS ANTHRACIS* SPORES

Michael Bartoszcze, Krzysztof Chomiczewski, and Lidia Mizak  
Military Institute of Hygiene and Epidemiology  
Lubelska 2, 24-100 Pulawy, Poland  
Tel/fax 48 818862822, e-mail obwwihe@man.pulawy.pl

### ABSTRACT

Fast and effective decontamination of areas contaminated with *Bacillus anthracis* spores after a successful bio-terrorist attack proves to be a challenging task. There exist a variety of disinfectants that can inactivate *Bacillus anthracis* spores; however, most of them have negative side effects, such as equipment corrosion and environmental toxicity. The investigation described here shows that Pera Safe™ has high sporocidal activity on *Bacillus anthracis* spores within 20 minutes, while being safe for use with minimal impact on the environment.

### INTRODUCTION

*Bacillus anthracis* is one of the main organisms that can be used in biological warfare or by bio-terrorists. This pathogen, in the form of spores, is characterized by considerable resistance to external factors, as is shown by its survival in the natural environment for years. *Bacillus anthracis* spores that are used as a biological weapon will result in environmental contamination that eventually must be decontaminated. Most of the disinfectants currently used to inactivate these spores are very corrosive, thus limiting the scope of their application. The purpose of this study was to estimate the sporocidal properties of a preparation of Pera Safe™ against *Bacillus anthracis* spores. A search of the literature showed no previous work with this material.

Report Documentation Page				Form Approved OMB No. 0704-0188	
Public reporting burden for the collection of information is estimated to average 1 hour per response, including the time for reviewing instructions, searching existing data sources, gathering and maintaining the data needed, and completing and reviewing the collection of information. Send comments regarding this burden estimate or any other aspect of this collection of information, including suggestions for reducing this burden, to Washington Headquarters Services, Directorate for Information Operations and Reports, 1215 Jefferson Davis Highway, Suite 1204, Arlington VA 22202-4302. Respondents should be aware that notwithstanding any other provision of law, no person shall be subject to a penalty for failing to comply with a collection of information if it does not display a currently valid OMB control number.					
1. REPORT DATE <b>00 JAN 2002</b>		2. REPORT TYPE <b>N/A</b>		3. DATES COVERED <b>-</b>	
4. TITLE AND SUBTITLE <b>Activity Of Pera Safe Against Bacillus Anthracis Spores</b>				5a. CONTRACT NUMBER	
				5b. GRANT NUMBER	
				5c. PROGRAM ELEMENT NUMBER	
6. AUTHOR(S)				5d. PROJECT NUMBER	
				5e. TASK NUMBER	
				5f. WORK UNIT NUMBER	
7. PERFORMING ORGANIZATION NAME(S) AND ADDRESS(ES) <b>Military Institute of Hygiene and Epidemiology Lubelska 2, 24-100 Pulawy, Poland</b>				8. PERFORMING ORGANIZATION REPORT NUMBER	
9. SPONSORING/MONITORING AGENCY NAME(S) AND ADDRESS(ES)				10. SPONSOR/MONITOR'S ACRONYM(S)	
				11. SPONSOR/MONITOR'S REPORT NUMBER(S)	
12. DISTRIBUTION/AVAILABILITY STATEMENT <b>Approved for public release, distribution unlimited</b>					
13. SUPPLEMENTARY NOTES <b>This article is from ADA409494 Proceedings of the 2001 ECBC Scientific Conference on Chemical and Biological Defense Research, 6-8 March , Marriott's Hunt Valley Inn, Hunt Valley, MD.</b>					
14. ABSTRACT					
15. SUBJECT TERMS					
16. SECURITY CLASSIFICATION OF:			17. LIMITATION OF ABSTRACT <b>UU</b>	18. NUMBER OF PAGES <b>4</b>	19a. NAME OF RESPONSIBLE PERSON
a. REPORT <b>unclassified</b>	b. ABSTRACT <b>unclassified</b>	c. THIS PAGE <b>unclassified</b>			

## MATERIALS AND METHODS

For this study, a suspension of *Bacillus anthracis* spores, strain Sterne 34 F2, was used. The Pera Safe™ preparation (series 2194) was manufactured by Antec International and was supplied by Naturan LTD, Warsaw, Poland.

**Spore suspension preparation:** 200 ml of media, consisting of brain- heart agar with the addition of yeast extract and MnSo<sub>4</sub> ( 0.1 g/l), was set in a level Roux bottle and was inoculated with 5 ml of a 24 hour broth bouillon culture of *Bacillus anthracis* (BA), strain Sterne 34 F2, and incubated for 72 hours at 37° C. The culture surface was then washed with Ringer's solution ( dilution 1:4), filtered through sterile gauze, and centrifuged five times at 5-7,000 revolutions per minute. After each centrifugation, the precipitate was washed with sterile PBS. The washed spores were suspended in sterile, distilled water until a titer of  $1.5 \times 10^7$  spores/ ml was obtained. The number of spores was determined in a Thoma chamber using phase-contrast microscopy.

**Preparation of the PeraSafe™ solution.** A solution of a 1.62% concentration was prepared in 500 ml flasks with sterile, distilled water at a of temperature 35° C. After dissolving, the solutions were left for 1 hour at room temperature. Immediately before testing, 1% bovine albumin (Serva 11925) was added to one of the flasks. This is a standard load of organic substances used in investigations of preparations for the disinfection of instruments.

**Determination of sporicidal activity.** Suspensions of BA spores with concentrations between  $1.5 \times 10^7$  spores/ml and  $1.7 \times 10^4$  spores/ml were used for testing the PeraSafe™ preparation. The suspensions were centrifuged, and then 10 ml of the PeraSafe™ solution was added to each of the precipitates, respectively.

Concurrently, the same investigation was performed using the sporicidal preparation with a 1% bovine albumin supplement. The samples were mixed thoroughly and left at room temperature for 20, 40, 80, 160 and 240 minutes, respectively. After the measured time elapsed, the suspension was centrifuged, washed, and the sporicidal activity was established by incubation of 100 ul of each of the samples on blood agar.

## RESULTS AND DISCUSSION

PeraSafe™ is a preparation in the form of powder containing sodium perborate, TAED, corrosion inhibitors, stabilizers, and dye. After dissolving the powder in water, supra-octane ions are liberated. This powder is used for general disinfection and shows bactericidal, fungicidal and virocidal properties.

The results are shown on Table 1. After 20 minutes, the solution of PeraSafe™ totally inactivates  $5 \times 10^4$  cfu/ml of *Bacillus anthracis* spores regardless of whether or not they are in the presence of bovine albumin. Higher concentrations of spores require a longer contact time, as is shown by the fact that a suspension of  $1.5 \times 10^6$  cfu/ml required 160 minutes for 100% inactivation. A concentration of  $1.5 \times 10^6$  spore/ml required at least 80 minutes of reaction time, irrespective of the presence of bovine albumin.

TABLE 1. Activity of Perasafe™ Against *Bacillus Anthracis* Spores.

PeraSafe™ contact time (min.)	Spore Concentration (cfu/ml)							
	1.5 x 10 <sup>7</sup>		1.5 x 10 <sup>6</sup>		1.5 x 10 <sup>5</sup>		1.5 x 10 <sup>4</sup>	
	Without albumin	With albumin	Without albumin	With albumin	Without albumin	With albumin	Without albumin	With albumin
20	3,0 x 10 <sup>2</sup>	5,4 x 10 <sup>2</sup>	2,7 x 10 <sup>1</sup>	4,0 x 10 <sup>1</sup>	0	2,0 x 10 <sup>1</sup>	0	0
40	1,1 x 10 <sup>2</sup>	2,1 x 10 <sup>2</sup>	1,2 x 10 <sup>1</sup>	2,0 x 10 <sup>1</sup>	0	0	0	0
80	2,7 x 10 <sup>1</sup>	5,0 x 10 <sup>1</sup>	0	0	0	0	0	0
160	0	0	0	0	0	0	0	0
240	0	0	0	0	0	0	0	0

The results suggest that the PeraSafe™ preparation can be considered as one of the best products acting against BA spores. Currently, the following preparations are typically used for decontamination: formaldehyde, glutaraldehyde, hydrogen peroxide, peroxyacetic acid, chloramine, chlorine water(1,2,3,4,5). However, chlorine solutions corrode metals, oxidize rubber, and are rapidly neutralized by organic substances. Formaldehyde and glutaraldehyde solutions are harmful for the skin and respiratory tract. Likewise, one needs long exposure times for effective decontamination with these preparations (up to a few hours), considerably longer than PeraSafe™.

Pera Safe acts rapidly and effectively, even in cases of considerable contamination with BA spores. This preparation is not only safe to use but also is characterized by a pleasant odor.

## CONCLUSIONS

PeraSafe™ inactivates *Bacillus anthracis* spores. The effectiveness of the sporicidal activity depends on the contact time and the concentration of spores. This preparation is safe for the environment.

## REFERENCES

1. Baldry M.G.C.: The bactericidal, fungicidal and sporicidal properties of hydrogen peroxide and peracetic acid. J.Appl. Bacteriol. 1983,54,417-23.
2. Dietz P., Böhm R.: Results of an experimental study on testing disinfectants with spores of *Bacillus anthracis*. Hyg. und Med. 1980, 5, 103-7.
3. Fair G.M., Morris J.C.: The dynamics of chlorination. J. New. Eng. Water Works Ass. 1947, 61, 285-301.
4. Gorman S.P., Scott E.M., Russell A.D.: A review – antimicrobial activity, uses and mechanisms of action of glutaraldehyde. J. Appl. Bacteriol. 1980, 48, 161-90.
5. Hussaini S.N., Ruby K.R.: Sporicidal activity of peracetic acid against *Bacillus anthracis* spores. Vet. Rec. 1976, 98, 257-9.